



BAUSCH + LOMB

Producto	Curva Base	Diam.	Material	Contenido de agua	Dk/t	Zona óptica	Rango de Graduaciones
Optima 38	8.70	14.00	Polymacon	38.6%	14	9-10mm	Esf. -0.25 a -.900 pasos de 0.25D
							Esf. -9.00 a -12.00 pasos de 0.50D
Soflens 59 caja c/6	8.60	14.20	Hilafilcon B	59%	16	9.0mm	Esf. -0.50 a -6.00 pasos de 0.25D Esf. -6.00 a -9.00 pasos de 0.50D Esf. '+0.50 a +6.00 pasos de 0.25D
Soflens 66 Toric caja c/6	8.50	14.50	Alphafilcon A	66%	16	8.0mm	Esf. +6.00 a-6.00 pasos 0.25D Esf. -6.00 a -9.00 pasos 0.50D Cil. -0.75, -1.25, -1.75, -2.25, -2.75 Ejes 10°, 20°, 90°, 160°, 170°, y 180°
							Neutro
Purevision 2 caja c/6	8.60	14.00	Balafilcon A	36%	130	9.0mm	Esf. +6.00 a -6.00 pasos 0.25D Esf. -6.00 a -9.00 pasos 0.50D incluye neutro
							Esf. +6.00 a -6.00 pasos 0.25D Esf. -6.00 a -9.00 pasos 0.50D Cil. -0.75, -1.25, -1.75, -2.25, -2.75 Ejes 10°, 20°, 90°, 160°, 170°, y 180°
Ultra caja c/6	8.50	14.20	Samfilcon A	46%	163	9.0mm	Esf. +6.00 a -6.00 pasos 0.25D Esf. -6.00 a -9.00 pasos 0.50D incluye neutro
Ultra Asigmatism caja c/6	8.60	14.50	Samfilcon A	46%	114	8.0mm	Esf. +6.00 a -6.00 pasos 0.25D Esf. -6.00 a -9.00 pasos 0.50D Cil. -0.75, -1.25, -1.75, -2.25, -2.75 Ejes 10° a 180° pasos de 10°
Ultra Multifocal caja c/6	8.50	14.20	Samfilcon A	46%	163	6.0mm	Esf. +6.00 a -10.00 pasos 0.25D Add. Low +0.75 a +15.0 Add High +1.75 a +2.50
Lunare caja c/2	8.70	14.10	Polymacon	38.6%	29	5.0mm	Neutro
Lunare caja c/1							Esf. 0.00 a -5.00 pasos 0.25D Esf. -5.00 a -6.00 pasos 0.50D

Soluciones para lentes de contacto	
Renu Fresh	60 ml.
	120 ml.
	355 ml.
	500 ml.
	355 ml. + 60ml.
Renu gotas rehumectantes	8ml.
Renu Lens case	
Simplus	105 ml.
Kit inicio Renu Fresh	
Biotrue	60 ml.
	120 ml.
	300 ml.
Renu Advanced	60ml.
	120ml.
	355ml.

La familia completa de lentes de contacto

ULTRA® contact lenses
with **MoistureSeal™** technology



16

HORAS DE COMODIDAD
con la tecnología MoistureSeal™

ULTRA® comodidad durante todo el día,
visión clara y consistente^{1,2}

CONSULE A SU PROFESIONAL DE LA SALUD VISUAL

Referencias: 1.. B+L unpublished data, 2013. 16-hour bilateral dispensing clinical evaluation of methacrylamide lenses compared to Air Optix Aqua, Acuvue Oasys and Biofinity. 2. Results of an online survey with patients who wore their lenses for 7+ days and on average, spend 8 or more hours a day on a digital device (n=45). Survey questions were top 3 box scores (% Strongly Agree, Agree, Slightly Agree) on a 6-point agreement scale, with a margin of error ±2.4%.

BAUSCH+LOMB

Reg. N o. 1348C20 16 SSA.